


The PHOTON 100: The First CMOS Active Pixel Sensor for X-ray Crystallography

Bruker AXS
Madison, WI



Welcome



Today's Topics

- Image Sensors
 - CMOS vs CCD
 - Scientific applications
- Introducing the new PHOTON 100
- D8 Crystallography Solutions
 - D8 QUEST
 - D8 VENTURE
 - PHOTON 100 data

Speakers



Michael Ruf
Global Product Manager, SC-XRD
Madison, WI USA





Roger Durst
Executive VP and CTO
Madison, WI USA

Image sensors

http://en.wikipedia.org/wiki/Image_sensor





- An image sensor is a device that converts an optical image into an electronic signal
- A modern sensor is typically a charge-coupled device (CCD) or a complementary metal-oxide-semiconductor (CMOS) active pixel sensor
- Both types of sensor accomplish the same task of capturing light and converting it into electrical signals
- Both types of sensor are silicon devices



CMOS versus CCD



- CMOS detector technology has begun to rapidly displace CCDs in most high-end applications
- CMOS technology provides many advantageous features and is quickly developing
- CCD technology has matured and has reached its performance limit "as good as it gets"



CMOS Sensor Technology in professional photography



- Google:
"Best professional digital SLR cameras review"
- All major manufacturers' flagship professional digital SLR cameras use large format CMOS sensors
 - **Nikon's** D3X uses Nikon's FX-format CMOS sensor
 - **Canon's** EOS-1D Mark III/IV uses Canon's full frame CMOS sensor
 - **Olympus' E-3** uses the 4/3 type Hi-Speed Live MOS sensor
 - **Sony** Alpha DSLR-A900
 - **Pentax** K-7



CMOS Sensor Technology in professional photography






- Current good quality sensors in digital cameras are photon noise limited (counting statistics)
- The foremost factor distinguishing high-end sensors for professional use from other image sensors is
 - The larger pixel size, which allows it to collect more light
 - Larger sensor size
- Larger pixels have higher signal-to-noise ratio at all levels, but especially at low signal levels
- The Nikon FX-format sensor has a size of 36 mm × 24 mm with a pixel size of 8.45 μm × 8.45 μm
- **Bruker's CMOS sensor is 11.5 times larger and its pixels have a 130 times larger area!**



CMOS Sensor Technology in scientific applications



- Astronomy - Stellar scintillation
 - the light of the star is refracted many times and in random directions when it hits a change in density in earth's atmosphere
 - Sensitive, low-noise, high-speed CMOS sensors are used to selectively assemble undistorted images




**Astronomy: Palomar 5m
"Lucky Imaging"**

CMOS Sensor Technology in scientific applications






- Terrain Mapping Cameras
 - CMOS Hyper Spectral mapping of lunar surface in visible and near infrared spectral regions
 - ...Very high dynamic range of the sensor helps in taking excellent features of dark polar sections of moon surface...

- Light microscopy
 - ...CMOS is unique in its ability to simultaneously offer ultra-low noise, extremely fast frame rates, wide dynamic range, high resolution and a large field of view...



Satellite Imaging: TMC




Microscopy: Fluorescently labeled fixed cell
(x60 oil objective)

New PHOTON 100 Detector for X-ray crystallography



- Large 100 cm² sensor for fast and efficient data collection
- Low-power consumption sensor for high-reliability
- Air-cooled for low maintenance
- High sensitivity
- No glass fiber taper for superior spatial accuracy
- Optimized pixel size and point spread function for superior signal
- Large pixel volumes for best light conversion

PHOTON 100 CMOS Sensor for X-ray crystallography




- First CMOS active pixel sensor for X-ray diffraction
- Large 10 x 10 cm² active area
- Seamless, no gaps or dead areas
- 96 micrometer pixels
- High sensitivity, near quantum limited
- Fast 0.5 sec readout for fast data collection
- Real-time hardware data correction
- No frame correlation needed
- No blooming or streaking for best data quality
- *Capable of continuous scan mode – shutter-less operation*

Introducing the PHOTON 100


- The first CMOS active pixel sensor (APS) detector for crystallography
- 100 x 100 mm² active area
- Taper-free: no demagnification
- Air-cooled
- Cu/Mo/Ag
- High sensitivity
- High reliability
 - 3-year detector warranty
- Fast readout
 - *Zero readout dead time mode*





- CCD: charge is shifted from one pixel to the next until it reaches the readout amplifier
 - Similar to a shift register


CMOS active pixel sensor




- CMOS: charge is readout from each pixel through a readout bus
 - Similar to a random access memory
- X-rays converted in phosphor screen (like CCD)

PHOTON 100 CMOS imager

- The world's largest production sensor for X-ray imaging
 - 100 x 100 mm²
 - 1024 x 1024 pixels
 - Exclusive to Bruker
- 4 times larger than the KAF 4320 CCD used in many detectors



PHOTON 100: taper-free design




- Large 100-cm² sensor for fast and efficient data collection
- Low-power-consumption sensor for high reliability: backed by our unique three-year detector warranty
- Air-cooled for low maintenance
- No glass fiber taper for superior spatial accuracy and high sensitivity
- Optimized pixel size and point-spread function for superior signal
- Large pixel volumes for best light conversion
- Fast readout for fast data collection
- No frame correlation needed
- No blooming or streaking
- Perfect match of the pixel size and the point spread of the high-resolution scintillator screens


- Like the APEX II, the PHOTON 100 does not require a fiber optic taper
 - Better sensitivity
 - Better spatial resolution
 - No distortion (better unit cells)
 - No zingers!

CMOS continuous scans: “shutter free”

- Unlike a CCD, a CMOS detector can be readout without closing the X-ray shutter
 - This eliminates the dead time between frames
 - Significantly decreases “wall time”, especially for short exposures
 - Also eliminates timing jitter (which for short frames can be the dominant error)
 - 0.5 sec jitter error 2%!




a) Conventional (CCD) readout with mechanical shutter



b) CMOS readout with electronic (rolling) shutter

How do continuous scans improve data?

- Hasegawa et al. (SPRING-8) studied the effect of continuous scans on data quality
 - 2-sec exposures, lysozyme
 - Reported 10-20% improvement in I/σ depending on exposure time
 - Equivalent to 20-40% increase in intensity (or exposure time)




How does the PHOTON 100 compare to conventional detectors?

- Detective Collection Efficiency (DCE) quantifies the contribution of detector size, sensitivity, and readout speed on performance
- DCE is proportional to the “wall clock” time to achieve equivalent data statistics

$$DCE = \frac{\Omega}{4\pi} DQE (1 - f_{dead})$$

size → $\frac{\Omega}{4\pi}$
 sensitivity → DQE
 speed → $(1 - f_{dead})$



PHOTON 100: A “smart” camera

- In CCD detectors corrections (dark, flood, bad pixel, and spatial) are applied *after* readout in software
 - Slows down data acquisition and processing
- In the PHOTON 100, all corrections are applied *inside the camera*, in real time in a high performance pipeline processor
 - Data output from camera is *fully corrected*
 - No added dead time
 - Flood field and dark files stored *in the camera* non-volatile memory
 - New floods and darks can easily be collected and stored (in the camera) by the user
 - Up to 125,000 frames of data can be stored in the camera memory (250 GB)

Other CMOS advantages: Resistance to blooming

- CCDs suffer “blooming” of strong signals
 - Typically seen as a degradation in the quality of low resolution data
- CMOS is immune to blooming
 - Handles strong signals better, delivers superior low resolution data





PHOTON 100 comparison to CCDs: The details

	PHOTON 100	Agilent ATLAS
Active area (mm ²)	100 x 100	100 x 100
Sensor size (mm ²)	100 x 100	50 x 50
Gain (electrons)	450	80*
Detective Collection Efficiency	0.1	0.05
Dynamic range	16,000	14,000*
Taper	No	Yes
Spatial Distortion	No	Yes
Cooling	Air	Water
Anti-blooming	Yes	No
Real-time data correction	Yes	No
Continuous scans (shutter-free)	Yes	No
Warranty (years)	3	1

*Measured



Other PHOTON advantages: reliability

- CMOS detectors exhibit higher reliability (that is, a longer Mean Time Between Failures) due to lower operating voltage
 - CCD: 30 V
 - CMOS: 3 V
- ***We are passing this reliability advantage on in the form of the first (and only) 3-year factory warranty for the PHOTON 100***

Summary

- The PHOTON 100 is the next generation of X-ray detector technology
- Large active area
 - No fiber taper
 - No dead areas
 - No zingers
- Air-cooled
- High sensitivity
- Intrinsic anti-blooming
- Smart camera: built in processor, memory
 - Real time corrections
- Shutter-free operation
- High reliability: 3-year warranty

D8 Crystallography Solutions



D8 QUEST



D8 VENTURE

D8 QUEST Configurations





- Single-source sealed tube and I μ S systems
- TRIUMPH optional
- FIXED-CHI
 - Manual track
 - Motorized track
- KAPPA
 - Motorized track
- PHOTON 100
- All air-cooled configurations
- Low-temperature device
 - KRYOFLEX II
 - Cryostream

D8 VENTURE Configurations




- Single-source sealed tube and I μ S systems
- Dual-source systems
 - Sealed tube – I μ S
 - TRIUMPH optional
 - I μ S - I μ S
- FIXED-CHI
 - Manual track
 - Motorized track
- KAPPA
 - Motorized track
- PHOTON 100
- All air-cooled configurations
- Low-temperature device
 - KRYOFLEX II
 - Cryostream

PHOTON 100 data




- Large field of view
 - 87° in 2Theta at 5 cm distance
 - 109° in 2Theta at 4 cm distance
- Homogenous background
- No gaps
- Low noise
 - Beamstop shadow and mount shadow are visible from a 10s exposure
- Large dynamic range
- Good spatial resolution

PHOTON 100 data
D8 VENTURE with I μ S™ Cu MX




- Diffraction pattern of a Li-Phthalocyanine complex
- Strong diffraction spots at 147° 2Theta (0.82 Å) 30s exposure
- Detector edge at 152° 2Theta (0.80 Å)
- No spots detected on a customer Mo-sealed tube system

PHOTON 100 data - Bispyrazolone D8 VENTURE with I μ S™ Cu MX



- Improved instrument design allows data to higher resolution
 - 0.81 Å D8 VENTURE
 - 0.84 Å X8 PROSPECTOR
- Good structure from small weakly diffracting organic sample
 - 0.06 mm x 0.11 mm x 0.14 mm
 - R1 = 4.45%

PHOTON 100 data – Ylid D8 QUEST with Mo TRIUMPH




- Ylid is the best understood structure at Bruker and invaluable for judging instrument quality
- Ylid meets manufacturing specifications
- PHOTON 100 handles strong and weak reflections well
- $R1 = 2.45\%$



PHOTON 100 data– Ruby D8 QUEST with Mo TRIUMPH



- D8 QUEST with PHOTON 100 delivers best unit cell data
- Ruby unit cell data is very accurate
 - Better than 1 in 1000 or 0.1%
- Handles strong diffraction well
 - $R(\text{int}) = 1.89\%$ $R(\text{sigma}) = 0.76\%$
 - $R1 = 2.05\%$ to 0.6 \AA



PHOTON 100 data – Vitamin C D8 VENTURE with I μ S™ Cu MX



- D8 VENTURE handles well diffracting samples better
- Longer exposure times are possible without saturating the detector
- Weaker signals can be detected more accurately
- Better data than X8 PROSPECTOR (APEX II CCD)
 - D8 VENTURE: R1 = 2.39%
 - X8 PROSPECTOR: R1 = 2.60%

Small Environmental Footprint Systems





- For the first time completely air-cooled systems are available
- No compromise on performance
- Preservation of natural resources
 - 1,700 m³ of water saved per year
 - 43,000 kWh of electrical energy per year
- Reduced carbon footprint
 - 26 metric tons of CO₂ per year
- Reduced noise
 - No recirculating chillers


Q & A



Any Questions?

Please **type** any questions you may have for our speakers in the **Q&A panel** and click **Send**.





Innovation with Integrity